


Anti-Rest/Nrsf Rabbit Polyclonal Antibody For If/Icc, Elisa - Q13127

Item Number: CM0024564



Introduction

Rabbit polyclonal antibody targeting human REST/NRSF (RE1-silencing transcription factor). Validated for IF/ICC and ELISA. Ideal for studying neuronal gene regulation, stem cell maintenance, and disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	REST/NRSF Rabbit pAb
Remarks/Alias	WT6; XBR; HGF5; NRSF; DFNA27; GINGF5; REST/NRSF
Species	Human
GeneID	5978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 760-1060 of human REST/NRSF (NP_001180437.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q13127
Protein Weight	122kDa
Shipping	Ice bag