


Anti-Rest/Nrsf Rabbit Polyclonal Antibody For Wb, Elisa - Q13127

Item Number: CM0024956



Introduction

Rabbit polyclonal antibody targeting REST/NRSF, a neuronal gene repressor, validated for WB and ELISA in human, mouse, and rat samples. Essential for chromatin modification and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	REST/NRSF Rabbit pAb
Remarks/Alias	WT6; XBR; HGF5; NRSF; DFNA27; GINGF5; REST/NRSF
Species	Human
Gene ID	5978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-130 of human REST/NRSF (NP_005603.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13127
Protein Weight	122kDa
Shipping	Ice bag