

Anti-Ret Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P07949

Item Number: CM0002652



Introduction

Rabbit monoclonal antibody against human Ret proto-oncogene. Validated for WB, IF/ICC, ELISA with human and mouse reactivity. Suitable for cell signaling, cancer, and neurodevelopmental research.

[Learn More](#)

Parameter	Value
Product Name	Ret Rabbit mAb
Remarks/Alias	PTC; MTC1; HSCR1; MEN2A; MEN2B; CDHF12; CDHR16; RET-ELE1; Ret
Species	Human
GeneID (Human)	5979
GeneID	5979
Immunogen	Synthetic peptide corresponding to amino acids 1000-1100 of human Ret (NP_066124.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07949
Protein Weight	124 kDa
Shipping	Ice bag