

Anti-Rev-Erb Beta/Nr1D2 Monoclonal Antibody For Wb, Elisa - Q14995

Item Number: CM0002713



Introduction

Rabbit monoclonal Rev-erb beta/NR1D2 antibody for WB and ELISA. Detects human, mouse, and rat Rev-erb beta, a circadian rhythm and metabolic regulator. Ideal for research on nuclear receptors and chronobiology.

[Learn More](#)

Parameter	Value
Product Name	Rev-erb beta/NR1D2 Rabbit mAb
Remarks/Alias	RVR; BD73; EAR-1R; REVERBB; REVERBbeta; Rev-erb beta/NR1D2
Species	Human
GeneID (Human)	9975
GeneID	9975
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human Rev-erb beta/NR1D2 (NP_005117.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14995
Protein Weight	65kDa
Shipping	Ice bag