

Anti-Rev-Erb Beta/Nr1D2 Polyclonal Antibody For Wb, Elisa - Q14995

Item Number: CM0020866



Introduction

High-quality rabbit polyclonal antibody against human Rev-erb beta/NR1D2, suitable for WB and ELISA with cross-reactivity to mouse and rat. Ideal for circadian rhythm and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	Rev-erb beta/NR1D2 Rabbit pAb
Remarks/Alias	RVR; BD73; EAR-1R; REVERBB; REVERBbeta; Rev-erb beta/NR1D2
Species	Human
GeneID	9975
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human Rev-erb beta/NR1D2 (NP_005117.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14995
Protein Weight	65kDa
Shipping	Ice bag