


Anti-Rev-Erb α /Nr1D1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P20393

Item Number: CM0002748



Introduction

High-quality rabbit monoclonal antibody against human Rev-erb α /NR1D1 (P20393). Validated for WB, IF/ICC, IF-P, and ELISA. Reacts with human, mouse, and rat. Ideal for circadian rhythm and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	Rev-Erb α /NR1D1 Rabbit mAb
Remarks/Alias	EAR1; hRev; THRA1; THRAL; ear-1; REVERBA; REVERBalpha; Rev-Erb α /NR1D1
Species	Human
GeneID (Human)	9572
GeneID	9572
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20393
Protein Weight	67kDa
Shipping	Ice bag