


Anti-Rexo2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y3B8

Item Number: CM0014891



Introduction

High-quality rabbit polyclonal antibody against human REXO2, validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying mitochondrial RNA degradation and transcription initiation.

[Learn More](#)

Parameter	Value
Product Name	REXO2 Rabbit pAb
Remarks/Alias	RFN; SFN; REX2; CGI-114; REXO2
Species	Human
GeneID (human)	25996
GeneID	25996
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-237 of human REXO2 (NP_056338.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3B8
Protein Weight	27kDa
Shipping	Ice bag