

Anti-Rfc1 Rabbit Polyclonal Antibody For Wb, If/Icc And Elisa - P35251

Item Number: CM0021212



Introduction

RFC1 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets replication factor C subunit 1 (128kDa), essential for PCNA loading in DNA replication - P35251.

[Learn More](#)

Parameter	Specification
Product Name	RFC1 Rabbit pAb
Remarks/Alias	A1; RFC; PO-GA; RECC1; CANVAS; MHCBBF; RFC140; RFC1
Species	Human
GeneID (Human)	5981
GeneID	5981
Immunogen	Synthetic peptide corresponding to amino acids 400-700 of human RFC1 (NP_001191676.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35251
Protein Weight	128kDa
Shipping	Ice bag