

Anti-Rfc3 Polyclonal Antibody For Wb And Elisa - P40938

Item Number: CM0030061



Introduction

Rabbit polyclonal antibody targeting human RFC3, validated for Western blot and ELISA applications. Detects RFC3 in human, mouse, and rat samples. Essential for DNA replication and PCNA loading research.

[Learn More](#)

Parameter	Value
Product Name	RFC3 Rabbit pAb
Remarks/Alias	RFC38, RFC3
Species	Human
GeneID (human)	5983
GeneID	5983
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human RFC3 (NP_002906.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40938
Protein Weight	41kDa
Shipping	Ice bag