

Anti-Rfc4 Rabbit Polyclonal Antibody For Wb - P35249

Item Number: CM0029890



Introduction

This anti-RFC4 rabbit polyclonal antibody is validated for WB, IF/ICC, IP, and ELISA in human and rat samples, making it ideal for DNA replication and PCNA loading studies.

[Learn More](#)

Parameter	Value
Product Name	RFC4 Rabbit pAb
Remarks/Alias	A1; RFC37; RFC4
Species	Human
GeneID (Human)	5984
GeneID	5984
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 204-363 of human RFC4 (NP_002907.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	P35249
Protein Weight	40kDa
Shipping	Ice bag