

Anti-Rfc5 Polyclonal Antibody For Wb, Elisa - P40937

Item Number: CM0028513



Introduction

High-quality rabbit polyclonal antibody targeting human RFC5, validated for WB and ELISA. Cross-reacts with mouse and rat RFC5. Suitable for DNA replication and PCNA loading studies.

[Learn More](#)

Parameter	Value
Product Name	RFC5 Rabbit pAb
Remarks/Alias	RFC36; RFC5
Species	Human
GeneID (human)	5985
GeneID	5985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 181-340 of human RFC5 (NP_031396.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40937
Protein Weight	38kDa
Shipping	Ice bag