

Anti-Rfk Rabbit Polyclonal Antibody For Wb And Elisa - Q969G6

Item Number: CM0029170



Introduction

High-quality RfK rabbit polyclonal antibody for Western blot and ELISA applications, suitable for human and mouse samples. Validated with recombinant protein immunogen (aa 1-155), detects riboflavin kinase, a key enzyme in FAD synthesis.

[Learn More](#)

Parameter	Value
Product Name	RfK Rabbit pAb
Remarks/Alias	RIFK; RfK
Species	Human
GeneID (human)	55312
GeneID	55312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-155 of human RfK (NP_060809.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q969G6
Protein Weight	18kDa
Shipping	Ice bag