

Anti-Rfpl2 Rabbit Polyclonal Antibody For Wb, Elisa - 075678

Item Number: CM0017257

Introduction

RFPL2 rabbit polyclonal antibody validated for WB and ELISA. Detects human and rat RFPL2 (RING finger protein 79, 42 kDa), a zinc-finger protein expressed in prostate and brain. Suitable for urogenital biology research. SWISS: 075678.

[Learn More](#)

Parameter	Value
Product Name	RFPL2 Rabbit pAb
Remarks/Alias	RNF79; RFPL2
Species	Human
Gene ID (Human)	10739
Gene ID	10739
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human RFPL2 (NP_006596.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	075678
Protein Weight	42kDa
Shipping	Ice bag