


Anti-Rfwd2/Cop1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nhy2

Item Number: CM0017275



Introduction

Rabbit polyclonal antibody against human RFWD2/COP1 (E3 ubiquitin-protein ligase), validated for WB and ELISA; targets amino acids 582-731. Suitable for ubiquitination and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RFWD2 Rabbit pAb
Remarks/Alias	FAP78; RFWD2; CFAP78; RNF200
Species	Human
GeneID	64326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 582-731 of human RFWD2 (NP_071902.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NHY2
Protein Weight	80kDa
Shipping	Ice bag