

Anti-Rfx2 Polyclonal Antibody For Wb, Elisa - P48378

Item Number: CM0019704

Introduction

RFX2 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human and mouse RFX2, a transcription factor key to spermatogenesis and ciliogenesis. Ideal for reproductive biology and gene regulation studies.

[Learn More](#)

Specification	Detail
Product Name	RFX2 Rabbit pAb
Remarks/Alias	RFX2
Species	Human
Gene ID	5990 (Human)
Immunogen	Recombinant fusion protein containing amino acids 295-450 of human RFX2 (NP_000626.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48378
Protein Weight	80kDa
Shipping	Ice bag