

Anti-Rfx3 Rabbit Polyclonal Antibody For Wb, Elisa - P48380

Item Number: CM0012927

Introduction

Rabbit anti-RFX3 polyclonal antibody validated for Western blotting and ELISA. Targets human transcription factor RFX3, a key regulator of ciliogenesis and pancreatic islet cell differentiation. Suitable for developmental and ciliopathy research.

[Learn More](#)

Parameter	Value
Product Name	RFX3 Rabbit pAb
Remarks/Alias	RFX3
Species	Human
GeneID (Human)	5991
GeneID	5991
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 280-430 of human RFX3 (NP_602304.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P48380
Protein Weight	84kDa
Shipping	Ice bag