

Anti-Rga1 (Gpa1) Polyclonal Antibody For Wb, Elisa - Q0DJ33

Item Number: CM0022062

Introduction

Rabbit polyclonal antibody targeting *Oryza sativa* RGA1 (GPA1), suitable for Western blot and ELISA. Validated in rice, this antibody aids in studying plant growth signaling and gibberellin/brassinosteroid pathways. - Q0DJ33

[Learn More](#)

Parameter	Value
Product Name	RGA1 Rabbit pAb
Remarks/Alias	D1; GA1; GPA1; RGA1; GP-alpha-1; OJ1005_D04.15
Species	<i>Oryza sativa</i>
GeneID	4338448
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-230 of <i>Oryza sativa</i> RGA1 (Q0DJ33)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q0DJ33
Protein Weight	44kDa
Shipping	Ice bag