


Anti-Rgap1 Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, Elisa - Q9H0H5

Item Number: CM0009817



Introduction

RGAP1 (Rac GTPase-activating protein 1) rabbit monoclonal antibody, KD validated, for WB, IHC-P, ELISA. Detects human, mouse, rat RGAP1 (71 kDa). Targets a key cytokinesis regulator involved in cell division and erythropoiesis.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] RGAP1 Rabbit mAb
Remarks/Alias	CYK4; CDAN3B; ID-GAP; HsCYK-4; MgcRacGAP; [KD Validated] RGAP1
Species	Human
Gene ID (Human)	29127
Gene ID	29127
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 533-632 of human RGAP1 (NP_037409.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0H5
Protein Weight	71kDa
Shipping	Ice bag