

Anti-Rgl2 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Gxw1

Item Number: CM0020341

Introduction

Rabbit polyclonal antibody against Arabidopsis thaliana DELLA protein RGL2 (RGA-like 2), validated for Western blot and ELISA applications. Ideal for gibberellin signaling, seed germination, and plant stress response research.

[Learn More](#)

Parameter	Value
Product Name	RGL2 Rabbit pAb
Remarks/Alias	RGA-like 2; T21P5.13; T21P5_13; RGL2
Species	Arabidopsis thaliana
GeneID (human)	5863
GeneID	821251
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 40-156 of arabidopsis thaliana RGL2 (NP_186995.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q8GXW1
Protein Weight	60kDa
Shipping	Ice bag