

Anti-Rgl3 Rabbit Polyclonal Antibody - Rgl3

Item Number: CM0024091

Introduction

Polyclonal antibody targeting Arabidopsis thaliana RGL3 protein, suitable for Western blot and ELISA. Raised in rabbit against recombinant fusion protein, with validated cross-reactivity to Arabidopsis thaliana.

[Learn More](#)

Parameter	Value
Product Name	RGL3 Rabbit pAb
Species	Arabidopsis thaliana
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of Arabidopsis thaliana RGL3 (NP_197251.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
Shipping	Ice bag