

Anti-Rgma Rabbit Monoclonal Antibody For Wb And Elisa - Q96B86

Item Number: CM0003065



Introduction

Rabbit monoclonal antibody against human Repulsive Guidance Molecule A (RGMA), validated for WB and ELISA with rat cross-reactivity. Useful for investigating nervous system development, neurite outgrowth inhibition, and BMP signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	RGMA Rabbit mAb
Remarks/Alias	RGM; RGMA
Species	Human
Gene ID	56963
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human RGMA (Q96B86).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q96B86
Protein Weight	49 kDa
Shipping	Ice bag