


Anti-Rgn Polyclonal Antibody For Wb, Elisa - Q15493

Item Number: CM0027370



Introduction

Rabbit polyclonal antibody targeting human regucalcin (RGN/SMP30). Validated for WB and ELISA, detects human, mouse, and rat. Ideal for calcium signaling and gluconolactonase research.

[Learn More](#)

Parameter	Value
Product Name	RGN Rabbit pAb
Remarks/Alias	RC; GNL; SMP30; HEL-S-41; RGN
Species	Human
GeneID (Human)	9104
GeneID	9104
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-299 of human RGN (NP_004674.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15493
Protein Weight	33kDa
Shipping	Ice bag