

Anti-Rgr Polyclonal Antibody For Wb, Elisa - P47804

Item Number: CM0027041

Introduction

Rabbit polyclonal antibody targeting human RGR (RPE-retinal G protein-coupled receptor), validated for Western blot and ELISA, with cross-reactivity to rat. Supports research on visual cycle, retinal pigment epithelium, and retinoid isomerization.

[Learn More](#)

Parameter	Value
Product Name	RGR Rabbit pAb
Remarks/Alias	RP44; RGR
Species	Human
GeneID (Human)	5995
GeneID	5995
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 156-295 of human RGR (NP_002912.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P47804
Protein Weight	32kDa
Shipping	Ice bag