

Anti-Rgs10 Rabbit Polyclonal Antibody For Wb And Elisa - O43665

Item Number: CM0026219



Introduction

Rabbit polyclonal antibody targeting human RGS10, suitable for Western blot and ELISA applications. Cross-reacts with mouse and rat; detects a ~21kDa protein. Ideal for GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	RGS10 Rabbit pAb
Remarks/Alias	RGS10
Species	Human
Gene ID	6001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human RGS10 (NP_001005339.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43665
Protein Weight	21kDa
Shipping	Ice bag