

Anti-Rgs13 Polyclonal Antibody For Wb, Elisa - 014921

Item Number: CM0017965

Introduction

RGS13 polyclonal antibody (rabbit, WB/ELISA validated) targets human RGS13 (19kDa), a regulator of G protein signaling. Cross-reacts with mouse and rat. Ideal for signal transduction research.

[Learn More](#)

Parameter	Specification
Product Name	RGS13 Rabbit pAb
Remarks/Alias	RGS13
Species	Human
GeneID (Human)	6003
GeneID	6003
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-159 of human RGS13 (NP_002918.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14921
Protein Weight	19kDa
Shipping	Ice bag