


## Anti-Rgs14 Polyclonal Antibody For Wb, Ihc-P, Elisa - O43566

Item Number: CM0030142



### Introduction

High-quality rabbit polyclonal antibody against human RGS14, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Targets a 61kDa protein; ideal for GPCR signaling and memory research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | RGS14 Rabbit pAb                                                                                                   |
| Remarks/Alias    | RGS14                                                                                                              |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 10636                                                                                                              |
| GeneID           | 10636                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 230-510 of human RGS14 (NP_006471.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O43566                                                                                                             |
| Protein Weight   | 61kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |