


Anti-Rgs2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P41220

Item Number: CM0020384



Introduction

High-quality rabbit polyclonal antibody against human RGS2, validated for WB, IHC-P, and ELISA. Ideal for G protein signaling and blood pressure regulation research.

[Learn More](#)

Parameter	Value
Product Name	RGS2 Rabbit pAb
Remarks/Alias	G058; RGS2
Species	Human
GeneID (Human)	5997
GeneID	5997
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human RGS2 (NP_002914.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P41220
Protein Weight	24kDa
Shipping	Ice bag