

Anti-Rgs20 Rabbit Polyclonal Antibody For Wb, Elisa - 076081

Item Number: CM0025850

Introduction

Anti-RGS20 rabbit polyclonal antibody raised against recombinant human RGS20 (aa 1-241). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Targets regulator of G protein signaling 20, a Gz-selective GTPase-activating protein involved in signal transduction inhibition.

[Learn More](#)

Parameter	Value
Product Name	RGS20 Rabbit pAb
Remarks/Alias	RGSZ1; ZGAP1; gz-GAP; g(z)GAP; RGS20
Species	Human
GeneID (Human)	8601
GeneID	8601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human RGS20 (NP_003693.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	076081
Protein Weight	44kDa
Shipping	Ice bag