

Anti-Rgs4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P49798

Item Number: CM0022549

Introduction

CamelBio offers a rabbit polyclonal antibody against human RGS4 (SWISS P49798), validated for WB, IF/ICC, and ELISA. Recognizes mouse and rat homologs. Suitable for investigating GPCR signaling and regulation.

[Learn More](#)

Parameter	Description
Product Name	RGS4 Rabbit pAb
Remarks/Alias	RGP4; SCZD9; RGS4
Species	Human
GeneID (Human)	5999
GeneID	5999
Immunogen	Synthetic peptide within amino acids 1-100 of human RGS4 (NP_005604.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49798
Protein Weight	23 kDa
Shipping	Ice bag