


Anti-Rgs5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O15539

Item Number: CM0030235



Introduction

Rabbit polyclonal antibody against human RGS5 (Regulator of G protein signaling 5). Validated for WB, IHC-P, ELISA; cross-reacts with mouse. Suitable for GPCR signaling research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	RGS5 Rabbit pAb
Remarks/Alias	MST092; MST106; MST129; MSTP032; MSTP092; MSTP106; MSTP129; RGS5
Species	Human
Gene ID (Human)	8490
Gene ID	8490
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	O15539
Protein Weight	21kDa
Shipping	Ice bag