

Anti-Rgs7 Polyclonal Antibody For Wb And Elisa - P49802

Item Number: CM0025518



Introduction

Rabbit polyclonal antibody targeting human RGS7 (regulator of G protein signaling 7). Validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for investigating GPCR-mediated signal transduction cascades.

[Learn More](#)

Parameter	Value
Product Name	RGS7 Rabbit pAb
Remarks/Alias	RGS7
Species	Human
GeneID (human)	6000
GeneID	6000
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49802
Protein Weight	58kDa
Shipping	Ice bag