


Anti-Rhbdd1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Teb9

Item Number: CM0014363



Introduction

RHBDD1 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat RHBDD1, a rhomboid protease involved in ERAD, apoptosis, and spermatogenesis. Ideal for research.

[Learn More](#)

Parameter	Value
Product Name	RHBDD1 Rabbit pAb
Remarks/Alias	RRP4; RHBDL4; RHBDD1
Species	Human
GeneID	84236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-315 of human RHBDD1 (NP_115652.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TEB9
Protein Weight	36kDa
Shipping	Ice bag