

Anti-Rhbdf1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Cc6

Item Number: CM0024121

Introduction

Anti-RHBDF1 polyclonal antibody (rabbit) validated for WB, IF/ICC, ELISA. Targets human inactive rhomboid protein 1 (iRhom1), cross-reacts with mouse and rat. Ideal for EGFR signaling and protein transport research.

[Learn More](#)

Parameter	Value
Product Name	RHBDF1 Rabbit pAb
Remarks/Alias	Dist1; hDist1; C16orf8; EGFR-R5; gene-89; gene-90; RHBDF1
Species	Human
GeneID (Human)	64285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-400 of human RHBDF1 (NP_071895.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96CC6
Protein Weight	97kDa
Shipping	Ice bag