


Anti-Rheb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15382

Item Number: CM0018212



Introduction

RheB Rabbit pAb detects human, mouse, and rat RheB protein (~20 kDa) in WB, IF/ICC, and ELISA applications. Ideal for mTOR signaling and cell growth research.

[Learn More](#)

Parameter	Value
Product Name	RheB Rabbit pAb
Remarks/Alias	RHEB2; RheB
Species	Human
GeneID (human)	6009
GeneID	6009
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15382
Protein Weight	20kDa
Shipping	Ice bag