


Anti-Rhoa Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P61586

Item Number: CM0013958



Introduction

High-quality rabbit polyclonal antibody targeting RhoA (ARHA) for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Essential for studying cytoskeleton organization, cell migration, and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	RhoA Rabbit pAb
Remarks/Alias	ARHA; ARH12; RHO12; EDFAOB; RHOH12; RhoA
Species	Human
Gene ID (Human)	387
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 94-193 of human RhoA (NP_001655.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61586
Protein Weight	22kDa
Shipping	Ice bag