


Anti-RhoB Polyclonal Antibody For Wb, Elisa - P62745

Item Number: CM0029618



Introduction

High-quality rabbit anti-RhoB polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat RhoB (~22 kDa). Suitable for apoptosis, cancer, and cell signaling research. - P62745

[Learn More](#)

Parameter	Value
Product Name	RhoB Rabbit pAb
Remarks/Alias	ARH6; ARHB; RHOH6; MST081; MSTP081; RhoB
Species	Human
Gene ID (Human)	388
Gene ID	388
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-190 of human RhoB (NP_004031.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62745
Protein Weight	22kDa
Shipping	Ice bag