


Anti-RhoB Rabbit Monoclonal Antibody For Wb And Elisa - P62745

Item Number: CM0007615



Introduction

Rabbit monoclonal antibody targeting human RhoB. Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for apoptosis, tumorigenesis, and cell trafficking research.

[Learn More](#)

Parameter	Value
Product Name	RhoB Rabbit mAb
Remarks/Alias	ARH6; ARHB; RHOH6; MST081; MSTP081; RhoB
Species	Human
Gene ID	388
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 70-150 of human RhoB (NP_004031.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62745
Protein Weight	22kDa
Shipping	Ice bag