

Anti-Rhobtb2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Byz6

Item Number: CM0019933

Introduction

High-quality rabbit polyclonal antibody against human RHOBTB2 (DBC2/p83). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for cell proliferation and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	RHOBTB2 Rabbit pAb
Remarks/Alias	p83, DBC2, DEE64, EIEE64, RHOBTB2
Species	Human
GeneID (Human)	23221
GeneID	23221
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-480 of human RHOBTB2 (NP_001153508.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYZ6
Protein Weight	83kDa
Shipping	Ice bag