

Anti-Rhoc Polyclonal Antibody For Wb, If/Icc, Elisa - P08134

Item Number: CM0020657



Introduction

Rabbit polyclonal antibody against human RhoC, validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat RhoC (~22 kDa). Ideal for signal transduction and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	RhoC Rabbit pAb
Remarks/Alias	H9; ARH9; ARHC; RHOH9; RhoC
Species	Human
GeneID (human)	389
GeneID	389
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08134
Protein Weight	22kDa
Shipping	Ice bag