

Anti-Rhoc Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - P08134

Item Number: CM0014664



Introduction

KO-validated rabbit polyclonal antibody targeting RhoC (P08134). Designed for WB and ELISA, with cross-reactivity to human and rat. Enables investigation of RhoC's role in signal transduction, cytoskeletal dynamics, and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RhoC Rabbit pAb
Remarks/Alias	H9; ARH9; ARHC; RHOH9; oC
Species	Human
GeneID (Human)	389
GeneID	389
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human RhoC (NP_786886.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P08134
Protein Weight	22kDa
Shipping	Ice bag