


Anti-Rhodopsin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08100

Item Number: CM0019980



Introduction

Rhodopsin rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Target the G protein-coupled photoreceptor essential for dim-light vision and phototransduction research.

[Learn More](#)

Parameter	Value
Product Name	Rhodopsin Rabbit pAb
Remarks/Alias	RP4; OPN2; CSNBAD1; Rhodopsin
Species	Human
GeneID (Human)	6010
GeneID	6010
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human Rhodopsin (NP_000530.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P08100
Protein Weight	39kDa
Shipping	Ice bag