


Anti-Rhogap (P190-B) Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q13017

Item Number: CM0007667



Introduction

Rabbit monoclonal antibody targeting RhoGAP (p190-B, ARHGAP5). Validated for WB, IP, ELISA in human, mouse, rat. Ideal for studies of Rho GTPase signaling, cell migration, and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	RhoGAP Rabbit mAb
Remarks/Alias	GF12; p190-B; RhoGAP5; p190BRhoGAP; RhoGAP
Species	Human
GeneID (Human)	394
GeneID	394
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1150-1300 of human RhoGAP (Q13017)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13017
Protein Weight	172 kDa
Shipping	Ice bag