


Anti-Rhogap Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13017

Item Number: CM0025734



Introduction

RhoGAP Rabbit Polyclonal Antibody (Q13017) validated for WB, IF/ICC, ELISA in human, mouse, rat. GTPase-activating protein for Rho family. Ideal for cell signaling and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	RhoGAP Rabbit pAb
Remarks/Alias	GFI2; p190-B; RhoGAP5; p190BRhoGAP; RhoGAP
Species	Human
GeneID (human)	394
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1292-1501 of human RhoGAP (NP_001164.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13017
Protein Weight	172kDa
Shipping	Ice bag