

Anti-Rhogdi Polyclonal Antibody For Wb, Elisa - P52565

Item Number: CM0014115



Introduction

KO Validated anti-RhoGDI polyclonal antibody for WB, ELISA. Reacts with human, mouse, rat. Ideal for studying Rho GTPase signaling and cell motility regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RhoGDI Rabbit pAb
Remarks/Alias	GDIA1; NPHS8; RHOGDI; RHOGDI-1; HEL-S-47e; DI
Species	Human
Gene ID (Human)	396
Gene ID	396
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-204 of human RhoGDI (NP_001172006.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52565
Protein Weight	23kDa
Shipping	Ice bag