

Anti-RhoH Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15669

Item Number: CM0024801

Introduction

Rabbit polyclonal antibody against human RHOH (RhoH, TTF), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Suitable for hematopoietic cell signaling, T-cell development, and mast cell function studies.

[Learn More](#)

Parameter	Value
Product Name	RHOH Rabbit pAb
Remarks/Alias	TTF; ARHH; RHOH
Species	Human
GeneID (Human)	399
GeneID	399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-191 of human RHOH (NP_004301.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15669
Protein Weight	21kDa
Shipping	Ice bag