


Anti-Rhot1/Miro1 Polyclonal Antibody For Wb, Elisa - Q8Ixi2

Item Number: CM0026150



Introduction

Rabbit polyclonal antibody against human RHOT1/MIRO1 (mitochondrial Rho GTPase 1), validated for WB and ELISA. Detects human, mouse, and rat. Recombinant immunogen (aa 310-590). Ideal for mitochondrial trafficking and fission research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] RHOT1/MIRO1 Rabbit pAb
Remarks/Alias	ARHT1; MIRO1; MIRO-1; RHOT1/MIRO1
Species	Human
Gene ID (Human)	55288
Gene ID	55288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-590 of human RHOT1/MIRO1 (NP_060777.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IXI2
Protein Weight	71kDa
Shipping	Ice bag