

Anti-Rhot2 Rabbit Polyclonal Antibody For Wb - Q8Ixi1

Item Number: CM0027877



Introduction

Rabbit polyclonal antibody targeting human RHOT2 (MIRO-2) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Studies mitochondrial trafficking and subcellular distribution. UniProt Q8IXI1.

[Learn More](#)

Parameter	Value
Product Name	RHOT2 Rabbit pAb
Remarks/Alias	RASL; ARHT2; MIRO2; MIRO-2; C16orf39; RHOT2
Species	Human
GeneID (Human)	89941
GeneID	89941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-550 of human RHOT2 (NP_620124.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IXI1
Protein Weight	68kDa
Shipping	Ice bag