


Anti-Rhou Polyclonal Antibody For Wb, Elisa - Q7L0Q8

Item Number: CM0020766



Introduction

Rabbit polyclonal antibody targeting human RHO (WRCH1, CDC42L1), validated for Western blot and ELISA applications. Involved in cell morphology, cytoskeletal organization, and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	RHO Rabbit pAb
Remarks/Alias	ARHU; G28K; WRCH1; hG28K; CDC42L1; RHO
Species	Human
GeneID (human)	58480
GeneID	58480
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7L0Q8
Protein Weight	28kDa
Shipping	Ice bag