

Anti-Ric8A Rabbit Polyclonal Antibody For Wb And Elisa - Q9Npq8

Item Number: CM0017595

Introduction

High-quality anti-RIC8A rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat RIC8A (synembryn-A), a chaperone and GEF for G alpha proteins. Ideal for cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	RIC8A Rabbit pAb
Remarks/Alias	RIC8; RIC8A
Species	Human
GeneID (Human)	60626
GeneID	60626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 362-531 of human RIC8A (NP_001273063.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPQ8
Protein Weight	60kDa
Shipping	Ice bag