

Anti-Rictor Polyclonal Antibody For Wb - Q6R327

Item Number: CM0027561



Introduction

RICTOR rabbit polyclonal antibody validated for WB and ELISA. Recognizes human RICTOR with rat cross-reactivity. This ~192 kDa protein is a key adaptor in mTORC2, mediating AKT phosphorylation and cell growth. - Q6R327

[Learn More](#)

Parameter	Value
Product Name	RICTOR Rabbit pAb
Remarks/Alias	PIA, AVO3, hAVO3, RICTOR
Species	Human
GeneID	253260
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1043-1142 of human RICTOR (NP_689969.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q6R327
Protein Weight	192kDa
Shipping	Ice bag